

1. The shield consists of two pieces:
The shield and the headband



2. Start by bending the top part of
the shield and crease to form a seam



3. Insert the tabs on the shield into the
slots.



4. Bend the tabs forward and crease
firmly to lock the tabs in place.




5. The main shield is now assembled



4. Starting from the inside of the mask,
insert the tabs through the inner slot



5. Feed tabs back through the outer
slots to secure them in place




4. Staple the connection if possible to
reinforce the joint



5. Push open and crease the vent holes
to the desired position.



6. Mask is complete. Disinfect prior
to use.



PLEASE DISINFECT APPROPRIATELY BEFORE USE